Instruction

EURIA-PP

Pancreatic Polypeptide radioimmunoassay
For in vitro diagnostic use only

Document No.E-23-0023-11

February, 2016

EURIA-PP

English:	page	2
Francais:	page	14
Espanol:	página	26
Deutsch:	Seite	38
Italiano:	pagina	50
Svenska	sida	62

REF RB 316

IVD

INTRODUCTION

Pancreatic polypeptide (PP) is synthesized as an amino-terminal moiety of a precursor peptide. PP isolated from pancreas has 36 amino acid residues with an amidated C-terminal tyrosine. PP is secreted by F-cells of the islets of Langerhans. PP is localized almost entirely in the pancreas although detectable levels throughout gastrointestinal tract have been reported. PP in human plasma is reported to exist in at least four different forms: PP 1-36, PP 3-36 and two unidentified forms.

PP is released into plasma during stimulation of meal. The physiological role of PP includes inhibition of stimulated gastric and pancreatic exocrine secretions and augmentation of insulin inhibited hepatic glucose production. These actions of PP are mediated by specific receptors. Receptor binding studies have shown that the intact C-terminal tyrosine amide is necessary for biological activity.

PRINCIPLE OF THE METHOD

The intended use of these reagents is for assay of PP in human serum. PP in serum is assayed without extraction by a competitive radioimmunoassay using a rabbit antiserum raised against bovine PP. PP in standards and samples compete with ¹²⁵I-labelled human PP in binding to the antibodies. ¹²⁵I-PP binds in a reverse proportion to the concentration of PP in standards and samples. Antibody-bound ¹²⁵I-PP is separated from the unbound fraction using the double antibody-polyethyleneglycol precipitation technique. The radioactivity of the precipitates is measured. Human, synthetic PP is used for standardization.

For professional use within a laboratory.

CLINICAL CONSIDERATIONS

The secretion of PP is stimulated by meal especially protein and fat. PP is also produced by endocrine active tumours in the pancreas and the gastrointestinal tract. These tumours often produce several peptide hormones in the combinations PP-VIP, PP-glucagon or PP-gastrin. Tumours with only PP-secretion have been reported. These tumours may occur at the WDHA or Verner-Morrison syndrome.

Elevated fasting levels of PP in serum are found at the occurrence of PP-producing tumours and endocrine tumours in the pancreas and in the gastrointestinal tract.

Normal level of PP in human serum: <100 pmol/L (fasting level obtained with this procedure).

PRECAUTIONS

For in vitro diagnostic use only.

As the regulations may vary from one country to another, it is essential that the person responsible for the laboratory are familiar with current local regulations, concerning all aspects of radioactive materials of the type and quantity used in this test.

This kit contains components of human origin. They have been tested by immunoassay for hepatitis B surface antigen, antibodies to HCV and for antibodies to HIV-1 and HIV-2 and found to be negative. Nevertheless, all recommended precautions for the handling of blood derivatives, should be observed.

This kit contains ^{125}I (half-life: 60 days), emitting ionizing X (28 keV) and γ (35.5 keV) radiations. Steps should be taken to ensure the proper handling of the radioactive material, according to local and/or national regulations. Only authorized personnel should have access to the reagents.

The following precautions should be observed when handling radioactive materials:

- Radioactive material should be stored in specially designated areas, not normally accessible to unauthorized personnel.
- Handling of radioactive material should be conducted in authorized areas only.
- Care should be exercised to prevent ingestion and contact with the skin and clothing. Do not pipette radioactive solutions by mouth.
- Drinking, eating or smoking should be prohibited where radioactive material is being used.
- Hands should be protected by gloves and washed after using radioactive materials.
- Work should be carried out on a surface covered by disposable absorbing material.
- Spills of radioactive material should be removed immediately, and all contaminated materials disposed as radioactive waste. Contaminated surfaces should be cleaned with a detergent.

The reagents in this kit contain sodium azide. Contact with copper or lead drain pipes may result in the cumulative formation of highly explosive azide deposits. On disposal of the reagents in the sewerage, always flush with copious amounts of water, which prevents metallic azide formation. Plumbing suspected of being contaminated with these explosive deposits should be rinsed thoroughly with 10% sodium hydroxide solution.

COMPOSITION OF THE REAGENT KIT

The reagents provided in each kit is sufficient for 100 tubes.

1. Anti-PP (Reagent A)

Rabbit antiserum raised against bovine PP. For 100 tubes. Lyophilized in 5.0 mL 0.5M phosphate buffer, pH 7.4, 2.5% human serum albumin and 0.5% NaN₃.

Reconstitution in 52 mL distilled water.

2. ¹²⁵I-PP (Reagent B)

Contains 28 KBq or 0.75 μCi of ¹²⁵I-hPP at the activity reference date. Produced by iodination of synthetic human PP. HPLC-purified, monoiodinated.

Specific activity: 1700-2100 μCi/nmol (62-77 MBq/nmol).

Lyophilized in 1.25 mL 0.5M phosphate buffer, pH 7.4, 2.5% human serum albumin, 0.5% NaN₃. Contains 0.12 mL normal rabbit serum.

Reconstitution in 12.5 mL distilled water.

3. Double antibody-PEG (Reagent C)

50 mL diluted goat anti-rabbit-Ig antiserum. Diluent: 0.05M phosphate buffer, pH 7.4, 0.25% human serum albumin and 0.05% NaN₃. Contains 7.5% polyethylene glycol 6000 (w/v).

4. Standard diluent (Reagent D)

10.0 mL PP-free human serum, lyophilized. Contains 500 KIU aprotinin (Trasylo[®] or equivalent)/mL. Reconstitution in 10.0 mL distilled water. For preparation of PP-working standards.

5. PP-standard, 2 000 pmol/L (8370 pg/mL) (Reagent E)

2.00 mL, 2 000 pmol/L synthetic human PP-standard. Lyophilized in 0.05M phosphate buffer, pH 7.4, 0.25% human serum albumin, 0.05% NaN₃. Reconstitution in 2.00 mL distilled water.

6. Assay buffer (Reagent F)

5.0 mL 0.05 M phosphate buffer, pH 7.4, 0.25% human serum albumin and 0.05% NaN₃.

To be used instead of antiserum in the non-specific binding test tubes.

7. Controls (Reagent G-H)

Lyophilized serum controls with low (G) and high (H) concentration of PP.

1.00 mL of each control after reconstitution. The PP concentrations of the controls are given on the label of the vials. Contains 0.05% NaN₃.

REAGENTS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

Distilled water

11-13x55 mm disposable tubes, polystyrene

Pipettes with disposable tips: 100 and 500 μl

Pipettes (glass): 1.00, 5.00 and 10.00 mL

Measuring cylinders: 25 mL and 50 mL

Vortex mixer

Centrifuge, refrigerated giving a minimum of 1700 x g

Gamma counter

SPECIMEN COLLECTION

Patients should be fasting 10 hours prior to sample collection.

Vein blood is collected in tubes without additives. The sample is allowed to clot. The serum is separated by centrifugation at +4° C. The serum should be frozen within 4 hours and stored at -18° C or lower until assayed. Repeated thawing and freezing should be avoided.

REAGENT PREPARATION AND STORAGE

Store all reagents at 2-8° C before reconstitution and use. The water used for reconstitution of lyophilized reagents should be distilled in an all-glass apparatus or be of corresponding purity. Dissolve the contents in a vial by gentle inversion and avoid foaming. The stability of the reagents is found on the labels of the vials. For lyophilized reagents the expiry dates are valid for the unreconstituted reagents. Reconstituted reagents are stable for 10 weeks (no longer than to the expiry date) stored correctly.

Reagent A: Anti-PP

Reconstitute with 52 mL distilled water.
Store at 2-8° C.

Reagent B: ¹²⁵I-PP

Reconstitute with 12.5 mL distilled water.
Store at -18° C or lower if reused.

Reagent C: Double antibody-PEG

Ready for use. Mix thoroughly before use.
Store at 2-8° C.

Reagent D: Standard diluent

Reconstitute with 10.0 mL distilled water.
Store at -18° C or lower if reused.

Reagent E: PP-standard, 2 000 pmol/L

Reconstitute with 2.00 mL distilled water.
Store at -18° C or lower if reused.
For preparation of PP-working standards, see radioimmunoassay procedure.

Reagent F: Assay buffer

Ready to use.
Store at 2-8° C.

Reagent G-H: Controls

Reconstitute with 1.00 mL distilled water. Store at -18° C or lower if reused.

RADIOIMMUNOASSAY PROCEDURE

Reconstitute the reagents as specified. Reagents should be brought to room temperature prior to use. Accuracy in all pipetting steps is essential. All tests (standards, controls and samples) should be performed in duplicate.

A complete assay includes:

Standards (St-tubes): 7 different concentrations; 0, 6.25, 12.5, 25.0, 50.0, 100 and 200 pmol/L.

Controls (C-tubes).

Samples (P-tubes).

Tubes for determination of the **non-specific binding (NSB-tubes).**

Tubes for determination of the **total radioactivity added (TOT-tubes).**

For an overview see table 1 on page 9.

PERFORMANCE

1. Reconstitute the reagents according to the instructions.
2. Prepare the PP-working standards by dilution of the PP-standard 2000 pmol/L (Reagent E) with the standard diluent (Reagent D) according to the following:

a/ 0.200 mL standard 2000 pmol/L	+ 1.800 mL diluent	= 200 pmol/L
b/ 1.00 mL standard 200 pmol/L	+ 1.00 mL diluent	= 100 pmol/L
c/ 1.00 mL standard 100 pmol/L	+ 1.00 mL diluent	= 50 pmol/L
d/ 1.00 mL standard 50 pmol/L	+ 1.00 mL diluent	= 25 pmol/L
e/ 1.00 mL standard 25 pmol/L	+ 1.00 mL diluent	= 12.5 pmol/L
f/ 1.00 mL standard 12.5 pmol/L	+ 1.00 mL diluent	= 6.25 pmol/L
g/ Standard diluent		= 0 pmol/L.
- Store the standard solutions at -18° C or lower if reused.
3. Pipette 100 µL of the standards (0-200 pmol/L), samples and controls in their respective tubes. Pipette 100 µL of the zero-standard in the NSB-tubes.
4. Pipette 500 µL anti-PP (Reagent A) to all tubes except the NSB- and TOT-tubes.
5. Add 500 µL assay buffer (Reagent F) to the NSB-tubes.
6. Vortex-mix and incubate for 20-24 hours at 2-8° C.
7. Pipette 100 µL ¹²⁵I-PP (Reagent B) to all tubes. The TOT-tubes are sealed and kept aside.
8. Vortex-mix and incubate for 20-24 hours at 2-8° C.
9. Pipette 500 µL double antibody-PEG (Reagent C) to all tubes except the TOT-tubes. Mix this reagent before pipetting.
10. Vortex-mix carefully and incubate for 30-60 minutes at 2-8° C.
11. Centrifuge the tubes for 15 minutes at +4° C (minimum 1700 x g).
12. Decant the supernatants immediately after centrifugation.
13. Count the radioactivity of the precipitates in a gamma counter (counting time: 2-4 minutes).

CALCULATION OF RESULTS

1. Subtract the average count rate (CPM) of the non-specific binding tubes from the count rates (CPM) of the replicates of standards, controls and samples.
2. A standard curve is generated by plotting the precipitated CPM, bound fraction in CPM or % B/TOT against the concentrations of the PP-standards. An example of a standard curve is given on page 10.
3. Interpolate the PP concentrations of the samples and controls from the generated standard curve.
4. The standard curve and the calculations of the concentrations in samples and controls can also be done by a computer method.

QUALITY CONTROL

In order to completely monitor the consistent performance of the radioimmunoassay there are some important factors which must be checked.

1. The found concentrations of the control sera

PP levels should be within the limits given on the labels of the vials.

2. Total counts

Counts obtained should approximate the expected CPM when adjusted for counter efficiency and radioactive decay. The content of ¹²⁵I-PP in this kit will give 10 500 CPM (-5, +20%) at activity reference date (counter efficiency = 80%).

3. Maximum binding (Bo/TOT)

Calculate for each assay the % bound radioactivity in the zero-standard: $\frac{B_0}{TOT} \times 100$

4. Non-specific binding (NSB/TOT)

Calculate for each assay the % non-specific binding: $\frac{NSB}{TOT} \times 100$

The % non-specific binding should be less than 7%.

5. Slope of standard curve

For example monitor the 80, 50 and 20% points of the standard line for run to run reproducibility.

ASSAY CHARACTERISTICS

Sensitivity

The lowest detectable concentration is 3 pmol/L. The figure corresponds to a decrease in binding of two x SD of the bound radioactivity in the zero-concentration standard.

Accuracy

A mean recovery of 104% (95-113%) was achieved when known amounts of hPP were added to human serum.

Precision

Intra assay variation

<u>Level</u>	<u>Coefficient of variation (%CV)</u>	<u>N</u>
28.8 pmol/L	2.6	10
108.5 pmol/L	1.8	10

Inter assay variation (total variation)

<u>Level</u>	<u>Coefficient of variation (%CV)</u>	<u>N</u>
38.8 pmol/L	2.0	10
99.3 pmol/L	3.5	10

Specificity

The following cross-reactions have been found

<u>Peptide</u>	<u>Cross-reaction</u>
Pancreatic polypeptide, human	100.0 %
Pancreatic polypeptide, bovine	120 %
Gastric inhibitory peptide, porcine	0.02 %
Cholecystokinin 39, porcine	0.02 %
Secretin, porcine	0.02 %
Gastrin 34, human	<0.01 %
Gastrin 17, human	<0.01 %
Glucagon, human, porcine	0.03 %
Insulin, porcine	<0.01 %
ACTH 1-39, porcine	<0.003%
Neuropeptide Y, human	<0.8 %
Peptide YY, human	<1.0 %

Interference


Samples displaying cloudiness, hemolysis, hyperlipemia or containing fibrin may give inaccurate results.

OUTLINE OF THE RIA PROCEDURE

Type of tubes	Tube no	Standard sample or control	Anti-PP (A)	Assay buffer (F)		¹²⁵ I-PP (B)		Double antibody-PEG (C)	
TOT	1-2	-	-	-	Vortex-mix and	100 µL	Vortex-mix and	-	Vortex-mix and
NSB	3-4	100 µL	-	500	incubate	100 µL	incubate	500 µL	and
Stand 0	5-6	100 µL	500 µL	-	for	100 µL	for	500 µL	incubate
Stand 6.25	7-8	100 µL	500 µL	-	20-24	100 µL	20-24	500 µL	for 30-60
Stand 12.5	9-10	100 µL	500 µL	-	hours at	100 µL	hours at	500 µL	min. at
Stand 25	11-12	100 µL	500 µL	-	2-8° C.	100 µL	2-8° C.	500 µL	2-8° C.
Stand 50	13-14	100 µL	500 µL	-		100 µL		500 µL	Centrifuge
Stand 100	15-16	100 µL	500 µL	-		100 µL		500 µL	15 min. at
Stand 200	17-18	100 µL	500 µL	-		100 µL		500 µL	1700 x g at
Control (G)	19-20	100 µL	500 µL	-		100 µL		500 µL	+4° C.
Control (H)	21-22	100 µL	500 µL	-		100 µL		500 µL	Decant and
Sample 1	23-24	100 µL	500 µL	-		100 µL		500 µL	count the
Sample 2	25-26	100 µL	500 µL	-		100 µL		500 µL	radioactivi-
etc.									ty of the
									precipitates.

Table 1

EXAMPLE OF PP STANDARD CURVE

REFERENCES / REFERENCES / REFERENCIAS / LITERATUR / BIBLIOGRAFIA / REFERENSER

1. Schwartz, T.W., Gingerich, R.L. and Tager, H.S.
Biosynthesis of pancreatic polypeptide: identification of precursor and cosynthesized product.
J Biol Chem 225:11494-11498, 1980.
2. Greider, M.H., Gersell, D.J. and Gingerich, R.L.
Ultrastructural localization of pancreatic polypeptide in the F cell of the dog pancreas.
J Histo Chem Society 26:1103-1108, 1978.
3. Gersell, R.J., Gingerich, R.L. and Greider, M.H.
Regional distribution and concentration of pancreatic polypeptide in human and canine pancreas.
Diabetes 28:11-15, 1979.
4. Chance, R.E., Moon, N.E. and Johnson, M.C.
Human pancreatic polypeptide (HPP) and bovine pancreatic polypeptide (BPP).
In B.M. Jaffe and H.R. Behlman (Eds).
Methods of hormone radioimmunoassay.
Academic Press, New York, 1979, 657-672.
5. Kimmel, J.R., Hayden, L.J. and Pollock, H.G.
Isolation and characterization of a new pancreatic polypeptide hormone.
J Biol Chem 250:9369-9376, 1975.
6. Adrian, R.E., Bloom, S.R., Bryant, M.G., Polak, J.M., Heitz, P.H. and Barnes, A.
Distribution and release of human pancreatic polypeptide.
Gut 17:940-944, 1976.
7. Hazelwood, R.L.
Synthesis, storage, secretion and significance of pancreatic polypeptide in vertebrates.
In S.J. Cooperstien and D. Watkins (Eds).
The islets of Langerhans, Academic Press, New York, 1981, p.p. 275-283.
8. Gingerich, R.L., Akpan, J.O., Leith, K.M. and Gilbert, W.R.
Patterns of immunoreactive pancreatic polypeptide in human plasma.
Regulatory Peptides 33:275-285, 1991.
9. Sun, Y.S., Brunicardi, F.C., Duck, P., Walfisch, S., Berlin, S.A., Chance, R.E.,
Gingerich, R.L., Elahi, D. and Andersen, D.K.
Reversal of abnormal glucose metabolism in chronic pancreatitis by administration of
pancreatic polypeptide.
Am J Surgery 151:130-140, 1986.
10. Seymour, N.E., Brunicardi, F.C., Chaiken, R.L., Lebovitz, H.E., Chance, R.E.,
Gingerich, R.L., Elahi, D. and Andersen, D.K.
Reversal of abnormal glucose production after pancreatic resection by pancreatic
polypeptide administration in man. Surgery 104:119-129, 1988.

SYMBOLS USED ON LABELS / SYMBOLES UTILISES SUR LES ETIQUETTES / SIMBOLOS UTILIZADOS EN LAS ETIQUETAS / ERLÄUTERUNG DER SYMBOLE. SIMBOLI USATI SULLE ETICHETTE / SYMBOLER SOM BRUKES PÅ ETIKETTER / SYMBOLER PÅ ETIKETTERNA.

	Batch code. Numéro de lot. Número de lote. Chargen-Nummer. Numero di lotto. Satsnummer.
	Catalogue number. Référence catalogue. Número de catálogo. Katalog-Nummer. Numero di catalogo. Katalognummer.
	Use-by date. Date de péremption. Fecha de caducidad. Verfallsdatum. Usare entro. Använd före.
	Temperature limit. Seuils de températures. Rango de temperature. Temperaturbereich. Limitazioni di temperatura. Förvaringstemperatur.
	Date of manufacture. Date de fabrication. Fecha de fabricacion. Datum der Herstellung. Data di produzione. Tillverkningsdatum.
	Contains radioactive substances. Contient des substances radioactives. Contiene sustancias radiactivas. Enthält radioaktive Stoffe. Contiene sostanze radioattive. Innehåller radioaktiva ämnen.
	Biological risk. Risque biologique. Riesgo biológico. Biologische Gefährdung. Rishio biologico. Biologisk risk.
	Consult instructions for use. Lire le mode d'emploi. Consulte las instrucciones de uso. Gebrauchsanweisung beachten. Leggere le istruzioni per l'uso. Läs instruktionsmanualen.
	In vitro diagnostic medical device. Dispositif médical de diagnostic in vitro. Producto sanitario para diagnóstico in vitro. In-vitro-Diagnostikum. Dispositivo medico-diagnostico in vitro. In vitro diagnostik medicinsk utrustning.
	Manufacturer. Fabricant. Fabricante. Hersteller. Produttore. Tillverkare.
	Contains sufficient for 100 tests. Contenu suffisant pour 100 tests. Contenido suficiente para 100 pruebas. Inhalt ausreichend für 100 Tests. Contenuto sufficiente per 100 test. Innehåller tillräckligt för 100 test.
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive. Conformément à la directive européenne 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro. La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnóstico in vitro. Konform mit Richtlinie 98/79/EG zu In-vitro-Diagnostika. Conformità alla direttiva 98/79/CE relativa ai dispositivi medico-diagnostici in vitro. Överensstämmer med direktiv 98/79/EG för medicintekniska produkter

REAG	A	Ab		Anti-PP. Anti-PP. Anti-PP. Anti PP. Anti PP. Anti-PP.
REAG	B	Ag	¹²⁵I	¹²⁵ I-PP. PP ¹²⁵ I. PP I- ¹²⁵ . ¹²⁵ I-PP. PP- ¹²⁵ I. ¹²⁵ I-PP.
REAG	C	DAB		Double antibody-PEG. PEG double anticorps. Doble anticuerpo con PEG. Doppel-Antikörper-PEG. Secondo anticorpo-PEG. Dobbelt antistoff, fast fase. Dubbel antikropp-PEG lösning.
REAG	D	DIL	CAL	Standard diluent. Diluant standard. Diluyente estándar. Standard Verdünner. Diluente dello standard. Standardspädningsmedel.
REAG	E	CAL	2000	PP standard 2000 pmol/L. Standard de PP, 2000 pmol/L. Estándar PP 2000 pmol/L. PP Standard 2000 pmol/L. PPstandard P 2000pmol/L. PP-standard 2000 pmol/L.
REAG	F	BUF	AS	Assay buffer. Tampon de dosage. Tampón de ensayo. Assaypuffer. Tampone. Spädningsbuffert.
REAG	G	CONTROL		Control, level 1 (low). Témoins, niveau 1 (bas). Control, nivel 1 (bajo). Kontrolle, Level 1 (niedrig). Controllo, livello 1 (normale). Kontroll, nivå 1 (låg).
REAG	H	CONTROL		Control, level 2 (high). Témoins, niveau 2 (élevé). Control, nivel 2 (alto). Kontrolle, Level 2 (hoch). Controllo, livello 2 (elevato). Kontroll, nivå 2 (hög).

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
EXAMPLE OF PP STANDARD CURVE

REFERENCES / REFERENCES / REFERENCIAS / LITERATUR / BIBLIOGRAFIA / REFERENSER

1. Schwartz, T.W., Gingerich, R.L. and Tager, H.S.
Biosynthesis of pancreatic polypeptide: identification of precursor and cosynthesized product.
J Biol Chem 225:11494-11498, 1980.
2. Greider, M.H., Gersell, D.J. and Gingerich, R.L.
Ultrastructural localization of pancreatic polypeptide in the F cell of the dog pancreas.
J Histo Chem Society 26:1103-1108, 1978.
3. Gersell, R.J., Gingerich, R.L. and Greider, M.H.
Regional distribution and concentration of pancreatic polypeptide in human and canine pancreas.
Diabetes 28:11-15, 1979.
4. Chance, R.E., Moon, N.E. and Johnson, M.C.
Human pancreatic polypeptide (HPP) and bovine pancreatic polypeptide (BPP).
In B.M. Jaffe and H.R. Behlman (Eds).
Methods of hormone radioimmunoassay.
Academic Press, New York, 1979, 657-672.
5. Kimmel, J.R., Hayden, L.J. and Pollock, H.G.
Isolation and characterization of a new pancreatic polypeptide hormone.
J Biol Chem 250:9369-9376, 1975.
6. Adrian, R.E., Bloom, S.R., Bryant, M.G., Polak, J.M., Heitz, P.H. and Barnes, A.
Distribution and release of human pancreatic polypeptide.
Gut 17:940-944, 1976.
7. Hazelwood, R.L.
Synthesis, storage, secretion and significance of pancreatic polypeptide in vertebrates.
In S.J. Cooperstien and D. Watkins (Eds).
The islets of Langerhans, Academic Press, New York, 1981, p.p. 275-283.
8. Gingerich, R.L., Akpan, J.O., Leith, K.M. and Gilbert, W.R.
Patterns of immunoreactive pancreatic polypeptide in human plasma.
Regulatory Peptides 33:275-285, 1991.
9. Sun, Y.S., Brunnicardi, F.C., Duck, P., Walfisch, S., Berlin, S.A., Chance, R.E.,
Gingerich, R.L., Elahi, D. and Andersen, D.K.
Reversal of abnormal glucose metabolism in chronic pancreatitis by administration of
pancreatic polypeptide.
Am J Surgery 151:130-140, 1986.
10. Seymour, N.E., Brunnicardi, F.C., Chaiken, R.L., Lebovitz, H.E., Chance, R.E.,
Gingerich, R.L., Elahi, D. and Andersen, D.K.
Reversal of abnormal glucose production after pancreatic resection by pancreatic
polypeptide administration in man. Surgery 104:119-129, 1988.

SYMBOLS USED ON LABELS / SYMBOLES UTILISES SUR LES ETIQUETTES / SIMBOLOS UTILIZADOS EN LAS ETIQUETAS / ERLÄUTERUNG DER SYMBOLE. SIMBOLI USATI SULLE ETICHETTE / SYMBOLER SOM BRUKES PÅ ETIKETTER / SYMBOLER PÅ ETIKETTERNA.

	Batch code. Numéro de lot. Número de lote. Chargen-Nummer. Numero di lotto. Satsnummer.
	Catalogue number. Référence catalogue. Número de catálogo. Katalog-Nummer. Numero di catalogo. Katalognummer.
	Use-by date. Date de péremption. Fecha de caducidad. Verfallsdatum. Usare entro. Använd före.
	Temperature limit. Seuils de températures. Rango de temperature. Temperaturbereich. Limitazioni di temperatura. Förvaringstemperatur.
	Date of manufacture. Date de fabrication. Fecha de fabricacion. Datum der Herstellung. Data di produzione. Tillverkningsdatum.
	Contains radioactive substances. Contient des substances radioactives. Contiene sustancias radiactivas. Enthält radioaktive Stoffe. Contiene sostanze radioattive. Innehåller radioaktiva ämnen.
	Biological risk. Risque biologique. Riesgo biológico. Biologische Gefährdung. Rishio biologico. Biologisk risk.
	Consult instructions for use. Lire le mode d'emploi. Consulte las instrucciones de uso. Gebrauchsanweisung beachten. Leggere le istruzioni per l'uso. Läs instruktionsmanualen.
	In vitro diagnostic medical device. Dispositif médical de diagnostic in vitro. Producto sanitario para diagnóstico in vitro. In-vitro-Diagnostikum. Dispositivo medico-diagnostico in vitro. In vitro diagnostik medicinsk utrustning.
	Manufacturer. Fabricant. Fabricante. Hersteller. Produttore. Tillverkare.
	Contains sufficient for 100 tests. Contenu suffisant pour 100 tests. Contenido suficiente para 100 pruebas. Inhalt ausreichend für 100 Tests. Contenuto sufficiente per 100 test. Innehåller tillräckligt för 100 test.
100	
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive. Conformément à la directive européenne 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro. La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnóstico in vitro. Konform mit Richtlinie 98/79/EG zu In-vitro-Diagnostika. Conformità alla direttiva 98/79/CE relativa ai dispositivi medico-diagnostici in vitro. Överensstämmer med direktiv 98/79/EG för medicintekniska produkter

REAG	A	Ab		Anti-PP. Anti-PP. Anti-PP. Anti PP. Anti PP. Anti-PP.
REAG	B	Ag	¹²⁵I	¹²⁵ I-PP. PP ¹²⁵ I. PP I- ¹²⁵ . ¹²⁵ I-PP. PP- ¹²⁵ I. ¹²⁵ I-PP.
REAG	C	DAB		Double antibody-PEG. PEG double anticorps. Doble anticuerpo con PEG. Doppel-Antikörper-PEG. Secondo anticorpo-PEG. Dobbelt antistoff, fast fase. Dubbel antikropp-PEG lösning.
REAG	D	DIL	CAL	Standard diluent. Diluant standard. Diluyente estándar. Standard Verdünner. Diluente dello standard. Standardspädningsmedel.
REAG	E	CAL	2000	PP standard 2000 pmol/L. Standard de PP, 2000 pmol/L. Estándar PP 2000 pmol/L. PP Standard 2000 pmol/L. PPstandard P 2000pmol/L. PP-standard 2000 pmol/L.
REAG	F	BUF	AS	Assay buffer. Tampon de dosage. Tampón de ensayo. Assaypuffer. Tampone. Spädningsbuffert.
REAG	G	CONTROL		Control, level 1 (low). Témoins, niveau 1 (bas). Control, nivel 1 (bajo). Kontrolle, Level 1 (niedrig). Controllo, livello 1 (normale). Kontroll, nivå 1 (låg).
REAG	H	CONTROL		Control, level 2 (high). Témoins, niveau 2 (élevé). Control, nivel 2 (alto). Kontrolle, Level 2 (hoch). Controllo, livello 2 (elevato). Kontroll, nivå 2 (hög).

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EURO DIAGNOSTICA AB