Instruction

## WIESLAB<sup>®</sup> ANCA screen kit

Enzyme immunoassay for detection of  
autoantibodies against Proteinase 3 (PR3-ANCA)  
and Myeloperoxidase (MPO-ANCA)

Break apart microtitration strips (12x8) 96 wells  
Store the kit at +2-8° C  
For in vitro diagnostic use only


Document No. E-23-0159-05  
May, 2015

### WIESLAB<sup>®</sup> ANCA screen kit

English: page .....	2
Français: page .....	12
Español: página .....	19
Deutsch: Seite .....	25
Italiano: pagina .....	32
Dansk: side.....	39
Svenska: sida.....	45

**REF** CP 111

**IVD**

## **INTENDED USE**

The Wieslab<sup>®</sup> ANCA screening test kit is an enzyme-linked immunosorbent assay (ELISA) for qualitative detection of IgG antibodies to proteinase 3 and myeloperoxidase (MPO) in human sera. The assay is used to detect antibodies in a single serum specimen. The results of the assay are to be used as an aid to the diagnosis of systemic vasculitis, especially Wegener's granulomatosis (WG) and microscopic polyangiitis (MP). The analysis should be performed by trained laboratory professionals. FOR IN VITRO DIAGNOSTIC USE.

**A positive result should always be confirmed by a semi-quantitative assay.**

## **Summary and explanation**

ANCAs (anti-neutrophil cytoplasmic antibodies) are a family of autoantibodies related to vasculitis and inflammatory disorders. Since 1985, when c-ANCA was shown to be related to Wegener's granulomatosis (1), interest in ANCAs has steadily increased, and today these antibodies are considered to be major diagnostic tools for the diagnosis and follow up of systemic vasculitis (13, 14).

The first method to detect ANCA was indirect immuno-fluorescence (IIF) performed on ethanol fixed granulocytes (2). This method yields two patterns, a cytoplasmic staining of the granulocyte denoting the presence of c-ANCAs, and a perinuclear staining denoting the presence of p-ANCAs. IIF was followed by ELISAs using the purified proteins (3,4).

The granulocyte is full of granules each with many different proteins. It was early shown that antibodies from systemic vasculitis patients bind to the alpha fraction containing the azurophil granules. The most important proteins were proteinase 3 (PR3)(5, 6) and myeloperoxidase (MPO)(7, 8). PR3 is a serine protease with a molecular weight of 29kD, and MPO is a dimer with a molecular weight of 140kD. Thus antibodies to proteinase 3 are termed PR3-ANCA, and antibodies to myeloperoxidase are termed MPO-ANCA.

Approximately 80-90% of WG patients manifest PR3-ANCA and 5-15% MPO-ANCA. One category of vasculitis is microscopic polyangiitis (MP). Most patients with active MP are characterised by positive ANCA test results, MPO-ANCA being more frequent than PR3-ANCA (9, 10, 11, 12).

## **Principle of the Wieslab<sup>®</sup> ANCA screen kit**

The wells of the microtiterstrip are coated with purified proteinase 3 and myeloperoxidase. During the first incubation, specific antibodies in diluted serum will bind to the antigen coating. The wells are then washed to remove unbound antibodies and other components. A conjugate of alkaline phosphatase labelled (goat) antibodies to human IgG binds to the antibodies in the wells in the second incubation. After a further washing step, detection of specific antibodies is obtained by incubation with substrate solution. The amount of bound antibodies correlates to the colour intensity and is measured as absorbance (optical density (OD)).

## **Warnings and precautions**

- The human serum components used in the preparation of the controls and calibrators in the kit have been tested for the presence of antibodies to human immunodeficiency virus 1 & 2 (HIV 1&2), hepatitis C (HCV) as well as hepatitis B surface antigen by FDA approved methods and found negative. Because no test methods can offer complete assurance that HIV, HCV, hepatitis B virus, or other infectious agents are absent, specimens and human-based reagents should be handled as if capable of transmitting infectious agents.
- The Centers for Disease Control and Prevention and National Institutes of Health recommended that potentially infectious agents be handled at the Biosafety Level 2.
- All solutions contain ProClin 300 as a preservative. Never pipette by mouth or allow reagents or patient sample to come into contact with skin. Reagents containing ProClin may be irritating. Avoid contact with skin and eyes. In case of contact, flush with plenty of water.

- Material safety data sheets for all hazardous components contained in this kit are available on request from Euro Diagnostica.

BUF	WASH	30X
DIL		
CONJ		

CONTROL	+
CONTROL	-
SUBS	pNPP

### Warning

Contains ProClin 300:

Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC no. 220-239-6] (3:1)

- H317: May cause an allergic skin reaction.  
 P264: Wash hands thoroughly after handling.  
 P280: Wear protective gloves/protective clothing/eye protection/face protection.  
 P302+352: IF ON SKIN: Wash with plenty of soap and water.  
 P333+313: If skin irritation or rash occurs: Get medical advice/attention.

### Specimen collection

The ANCA screen assay is for use with serum. Handle as if capable of transmitting infectious agents. Avoid using sera which are icteric, lipemic and hemolyzed.

Heat-inactivated sera can yield unspecific reactivities and should not be used.

Store serum between 2-8° C if testing will take place within five days. If specimens are to be kept for longer periods, store at -20° C or colder.

Do not use a frost-free freezer because it may allow the specimens to go through freeze-thaw cycles and degrade antibody. Samples that are improperly stored or are subjected to multiple freeze-thaw cycles may yield spurious results.

The NCCLS provides recommendations for storing blood specimens, (Approved Standard-Procedures for the Handling and Processing of Blood Specimens, H18A, 1990).

### Kit components and storage of reagents

- One frame with breakapart strips (12x8) wells), coated with purified proteinase 3 (6 strips red coloured) and myeloperoxidase (6 strips green coloured), one lid sealed in a foil pack with a dry pack.
- 0.75 mL negative control (NC) containing human serum in diluent (green colour).
- 0.75 mL positive control (PC) containing human serum in diluent (red colour).
- 13 mL conjugate containing alkaline phosphatase-labelled (Goat) antibodies to human IgG (blue colour).
- 32 mL Diluent (Dil) containing PBS (reddish colour).
- 13 Substrate pNPP.
- 30 mL Wash solution 30x concentrated.

All reagents in the kit are ready for use except wash solution and should be stored at + 2-8° C. Remove only the number of wells needed for testing, resealing the aluminium package carefully.

### Materials or equipment required but not provided

- Microplate reader with filter 405 nm.
- Precision pipettes with disposable tips.
- Washer for strips, absorbent tissue, tubes and a timer.

**PROCEDURE**

All solutions should be used at room temperature.

Remove only the number of wells needed for testing, resealing the aluminium package carefully.

**If the assay is performed without a shaker, incubation times should be adjusted as follows: serum incubation 10 minutes, no change; conjugate incubation 20-30 minutes; substrate incubation 20-30 minutes.**

**Preparation of washing solution**

In case salt crystals are observed in the vial with concentrated wash solution, place the vial at 37°C water bath until the crystals have dissolved before dilution of wash solution. Dilute 10 mL of the 30x concentrated wash solution in 290 mL distilled water. When stored at 2-8° C, the diluted wash solution is stable until the date of expiration of the kit.

**Dilution of serum and incubation**

1. Pipet 75 µL diluent (Dil) into each well.
2. Pipet 25 µL negative control (NC) into wells A1 and A2.
3. Pipet 25 µL positive control (PC) into wells B1 and B2.
4. Pipet 25 µL of patient serum no.1 into wells C1 and C2, serum no.2 into wells D1 and D2, serum no. 3 into wells E1 and E2 etc. according to the scheme below. Put on the lid and incubate for 10 minutes on a shaker.

	1	2
	<b>Column 1 (PR3-ag) Red</b>	<b>Column 2 (MPO-ag) Green</b>
<b>A</b>	75 µL Dil + 25 µL NC	75 µL Dil + 25 µL NC
<b>B</b>	75 µL Dil + 25 µL PC	75 µL Dil + 25 µL PC
<b>C</b>	75 µL Dil + 25 µL pat 1	75 µL Dil + 25 µL pat 1
<b>D</b>	75 µL Dil + 25 µL pat 2	75 µL Dil + 25 µL pat 2
<b>E</b>	75 µL Dil + 25 µL pat 3	75 µL Dil + 25 µL pat 3
<b>F</b>	75 µL Dil + 25 µL pat 4	75 µL Dil + 25 µL pat 4
<b>G</b>	75 µL Dil + 25 µL pat 5	75 µL Dil + 25 µL pat 5
<b>H</b>	75 µL Dil + 25 µL pat 6	75 µL Dil + 25 µL pat 6

**After serum incubation**

Wash 4 times with washing solution, filling and emptying the wells each time; after the last wash, empty the wells by tapping the strip on an absorbent tissue.

**Adding conjugate**

Add 100 µL conjugate to each well. Put on the lid and incubate for 10 minutes on a shaker (without a shaker, incubate 20-30 minutes).

**After conjugate incubation**

Wash 4 times with washing solution, filling and emptying the wells each time; after the last wash, empty the wells by tapping the strip on an absorbent tissue.

**Adding substrate**

Add 100 µL substrate pNPP to each well, incubate for 10-20 minutes (without a shaker, incubate 20-30 minutes). Read the absorbance at 405 nm on a microplate reader at 10-20 minutes after substrate addition and 20-30 minutes after substrate addition without a shaker.

**Calculations**

An optical density (OD) ratio for each patient sample is calculated as follows:

$$\text{OD ratio} = \frac{\text{OD patient sample}}{\text{OD of negative control}}$$

The patient sample is negative if the OD ratio is < 3.0. Equivocal if the OD ratio is 3.0-4.0 and positive if the OD ratio is > 4.0.

### Quality Control

The OD for the negative control should be  $< 0.3$

The OD for the positive control should be  $> 0.9$

If any of the values are not within their respective ranges, the test should be considered invalid and the test should be repeated.

The patient sample may contain unspecific reactivities if it is positive on both antigens and must be reassayed in a quantitative assay. Sometimes unspecific reactivities will disappear with a higher dilution.

The negative and positive controls are intended to monitor for substantial reagent failure. The positive control will not ensure precision at the assay cut-off. It is recommended that an additional control be assayed that is close to the cut-off. Additional controls may be tested according to guidelines or requirements of local, state, and/or federal regulations or accrediting organisations. Refer to NCCLS C24-A for guidance on appropriate QC practices.

### Interpretation of results

A sample with a OD ratio of :

$< 3.0$  = **Negative**       $3.0 - 4.0$  = **Equivocal**; Retest, if still equivocal retest by an alternative method  
 $> 4.0$  = **Positive**

Owing to the high sensitivity of the screening kit, in a few percent of cases samples with an OD ratio  $> 3$  may be found to be negative in a quantitative assay.

**An equivocal or a positive test result should always be confirmed by a quantitative assay.**

### Limitations

The individual patient's OD ratio can not be used as a measure of disease severity, as antibodies from different patients may differ from each other in affinity. Thus, it is difficult to obtain an absolute standardization of results.

The test should not be relied upon as the sole basis of decisions on clinical therapy, but should be used in combination with clinical symptoms and the results of other available tests.

Sera from patients with other autoimmune diseases and from normal individuals may contain potentially cross-reactive autoantibodies. Some individuals may be positive, with little or no evidence of clinical disease. On the other hand, some patients with active disease may have undetectable levels of these antibodies.

Immunosuppressive therapy should not be started on basis of a positive ANCA result. Initiation or changes in treatment should not be based on changes in ANCA concentration alone, but rather on careful clinical observation.

### Expected results

MPO and PR3-ANCA are rarely found in normal healthy individuals. The ANCA screening kit was tested with 131 normal sera. 131 were found to be negative. One new patient with PR3-ANCA is expected per 100000 individuals per year. Around 10% of patients with WG are negative in both IIF and ELISA. MPO-ANCA are found in approximately one half of the sera of MP (microscopic polyangiitis) patients. The MPO-ANCA was tested with 42 WG sera, 4 were found to be positive. The MPO-ANCA was tested with 43 MP sera, 20 were found to be positive. (Table 1)

The PR3-ANCA was tested with 42 WG sera, 39 were found to be positive.

PR3-ANCA are found in approximately one half of the sera of MP (microscopic polyangiitis) patients.

The PR3-ANCA was tested with 43 MP sera, 21 were found to be positive. (Table 1)

## Performance characteristics

**Table 1. Clinical sensitivity and specificity.** A total of 288 frozen retrospective sera with clinical characterisation were assayed. The following table summarises the results

Control and Disease groups	Total	Negative < 3		Equivocal 3-4		Positive >4	
		PR3	MPO	PR3	MPO	PR3	MPO
Blood donors: (NS)	131	131	127	0	4	0	0
WG:	42	3	37	0	1	39	4
MP:	43	20	23	2	0	21	20
SLE:	31	31	24	0	2	0	*5
RA:	41	41	40	0	1	0	0

WG = Wegener's granulomatosis, MP = microscopic polyangiitis RA = rheumatoid arthritis  
SLE = systemic lupus erythematosus GBM = glomerular basement membrane

\* All 5 sera was positive in quantitative MPO-ANCA test.

### Clinical sensitivity (Equivocal samples not included in the calculations)

#### PR3-ANCA

WG = 39/42 = 92.9 % 95% CI = 80.5-98.5% MP=21/41=51.2 % 95% CI = 35.1-67.1%

#### MPO-ANCA

WG = 4/41 = 9.7 % 95% CI = 2.7-23.1% MP=20/43=46.5 % 95% CI = 31.2-62.3%

### Clinical specificity (Equivocal samples not included in the calculations)

#### PR3-ANCA

SLE = 31/31 = 100 % 95% CI = 88.8-100%  
RA=41/41=100 % 95% CI = 91.4-100% NS = 131/131 = 100 % 95% CI = 97.2-100%

#### MPO-ANCA

SLE =24/29 =82.8 % 95% CI = 64.2-94.2%  
RA=40/40=100 % 95% CI = 91.2-100% NS = 127/127 = 100 % 95% CI = 97.1-100%

**Table 2. Relative sensitivity and specificity** of the Wieslab® PR3-ANCA screen kit compared to a semi-quantitative ELISA. A total of 216 frozen retrospective sera were assayed. The table summarises the results.

Wieslab® ANCA screen							
Semi-quantitative ELISA		Negative < 3		Equivocal 3-4		Positive >4	
		PR3	MPO	PR3	MPO	PR3	MPO
PR3-ANCA	Positive	1	1	0	0	59	23
MPO-ANCA	Negative	152	182	2	4	0	0
	Equivocal	1	4	0	1	1	1
	Total	154	187	2	5	60	24

### Sera falling in the equivocal range are not included in the calculations:

Relative sensitivity PR3-ANCA = 59/60 = 98.3 % 95% CI = 91.1-100%  
Relative sensitivity MPO-ANCA = 23/24 = 95.8 % 95% CI = 78.9-99.9%

Relative specificity PR3-ANCA = 152/152 = 100 % 95% CI = 97.6-100%  
Relative specificity MPO-ANCA = 182/182 = 100 % 95% CI = 98-100%

The 95% confidence interval (CI) was calculated using the exact method.

**Table 3. Batch to batch variation** was determined by testing three different samples in 4 different batches.

Sample		OD ratio		SD		CV %	
PR3	MPO	PR3	MPO	PR3	MPO	PR3	MPO
2	3	34.5	16.3	1.0	3.9	3	24
5	6	19.8	27.8	2.2	2.1	11	7
8	9	33.8	23.0	1.5	2.9	4	13

**Table 4. Inter-assay precision** was determined by testing one sample. Results were obtained for six different runs.

Sample		OD ratio		SD		CV %	
PR3	MPO	PR3	MPO	PR3	MPO	PR3	MPO
PK	PK	27.3	17.3	2.9	3.8	11	21
5	6	12.1	16.5	1.2	0.8	10	5

**Table 5. Intra-assay precision** was determined by testing one sample in 22 wells.

Sample		Mean OD		SD		CV %	
PR3	MPO	PR3	MPO	PR3	MPO	PR3	MPO
PK	PK	1.3	1.9	0.07	0.06	6	3
5	6	1.3	1.5	0.06	0.07	5	5

**TROUBLESHOOTING**

<b>Problem</b>	<b>Possible causes</b>	<b>Solution</b>
Control values out of range.	<ol style="list-style-type: none"> <li>1. Incorrect temperature, timing or pipetting; reagents not mixed.</li> <li>2. Cross contamination of controls.</li> <li>3. Improper dilution.</li> <li>4. Optical pathway not clean.</li> </ol>	<ol style="list-style-type: none"> <li>1. Check that the time and temperature was correct. See "Poor precision" below. Repeat test.</li> <li>2. Pipette carefully.</li> <li>3. Repeat test.</li> <li>4. Check for dirt or air bubbles in the wells. Wipe bottom and reread.</li> </ol>
All test results negative.	<ol style="list-style-type: none"> <li>1. One or more reagents not added, or added in wrong sequence.</li> <li>2. Antigen coated plate inactive.</li> </ol>	<ol style="list-style-type: none"> <li>1. Recheck procedure. Check for unused reagents. Repeat test.</li> <li>2. Check for obvious moisture in unused wells. Wipe bottom and reread.</li> </ol>
All test results yellow.	<ol style="list-style-type: none"> <li>1. Contaminated buffers or reagents.</li> <li>2. Washing solution contaminated.</li> <li>3. Improper dilution of serum.</li> </ol>	<ol style="list-style-type: none"> <li>1. Check all solutions for turbidity.</li> <li>2. Use clean container. Check quality of water used to prepare solution.</li> <li>3. Repeat test.</li> </ol>
Poor precision.	<ol style="list-style-type: none"> <li>1. Pipette delivery CV greater than 5%.</li> <li>2. Serum or reagents not mixed sufficiently or not equilibrated to room temperature.</li> <li>3. Reagent addition taking too long; inconsistency in timing intervals.</li> <li>4. Optical pathway not clean.</li> <li>5. Washing not consistent; trapped bubbles; washing solution left in the wells.</li> <li>6. Improper pipetting.</li> </ol>	<ol style="list-style-type: none"> <li>1. Check calibration of pipette. Use reproducible technique.</li> <li>2. Mix all reagents gently but thoroughly and equilibrate to room temperature.</li> <li>3. Develop consistent uniform technique and use multi-tip device or auto-dispenser to decrease time.</li> <li>4. Check for air bubbles in the wells. Wipe bottom and reread.</li> <li>5. Check that all wells are filled and aspirated uniformly. Dispense liquid above level of reagent in well. After the last wash, empty the wells by tapping the strip on an absorbent tissue.</li> <li>6. Avoid air bubbles in pipette tips.</li> </ol>

**REFERENCES. RÉFÉRENCES. REFERENCIAS. LITERATUR. BIBLIOGRAFIA. LITTERATUR. REFERANSER. REFERENSER**

- 1. van der Woude FJ, Rasmussen N, Lobatto S, Wiik A, Permin H, Van Es L et al.** Autoantibodies against neutrophils and monocytes: Tool for diagnosis and marker of disease activity in Wegener's granulomatosis. *Lancet* 1985, i 425-429.
- 2. Wiik A, Rasmussen N, Wieslander J.** Methods to detect autoantibodies to neutrophilic granulocytes, in Van Venrooij WJ, Maini RN (eds): *Manual of biological markers of disease*. Kluwer Academic Publishers 1993, A9, 1-14.
- 3. Rasmussen N, Sjölin C, Isaksson B, Bygren P, Wieslander J.** An ELISA for the detection of anti-neutrophil cytoplasm antibodies (ANCA). *J Immunol Methods* 1990, 127, 139-145.
- 4. Segelmark M, et al.** How and why should we detect ANCA ? *Clin Exp Rheumatol* 2000, 18, 629-635.
- 5. Jenne DE, Tschopp J, Lüdemann J, Utecht B, Gross WL.** Wegeners autoantigen decoded. *Nature* 1990, 346, 520.
- 6. Wieslander J, Wiik A.** ANCA antigens: Proteinase 3, in van Venrooij W, Maini R (eds): *Manual of Biological Markers of Disease*. Kluwer Academic Publishers 1994, 1-9.
- 7. Falk RJ, Jennette JC.** Anti-neutrophil cytoplasmic autoantibodies with specificity for myeloperoxidase in patients with systemic vasculitis and idiopathic necrotizing and crescentic glomerulonephritis. *N Engl J Med* 1988, 318, 1651-1657.
- 8. Wieslander J, Wiik A.** ANCA antigens: Myeloperoxidase, in van Venrooij W, Maini R (eds): *Manual of Biological Markers of Disease*. Kluwer Academic Publishers 1994, 1-9.
- 9. Kamech L, Harper L, Savage C.** ANCA - Positive Vasculitis. *I. Am Soc Nephrol* 2002 43 1953 - 1960.
- 10. Segelmark M, Elzouki AN, Wieslander J, Eriksson S.** Heterozygosity for the alpha 1-antitrypsin PiZ gene affects the outcome of PR3-ANCA positive vasculitis. *Kidney Int* 1995, 48, 844-850.
- 11. Boomsma MM, Stegeman CA, van der Leij MJ, et al.** Prediction of relapses in Wegener's granulomatosis by measurement of ANCA levels: a prospective study. *Arthritis Rheum*, 2000, 43 2025 - 2033.
- 12. Savige J, Davies D, Falk RJ, Jennette JC, and Wiik A.** ANCA and associated diseases: A review of the clinical and laboratory features. *Kidney Int* 2000, 57, 846-862.
- 13. Savige J, Gills D, Benson E et al.** International Consensus Statement on Testing and Reporting of ANCA. *Am J Clin Pathol* 1999, 111, 507-513.
- 14. Won H, et al.** Serial ANCA titers. Useful tool for the prevention of relapses in ANCA associated vasculitis. *Kidney Int*, 2003, 63, 1079-1085.

**Explanation of symbols. L'explication de symboles. La explicación de símbolos. Erklärung der Symbole. La spiegazione di simboli. Forklaring til symboler. Förklaringar till symboler.**

	Batch code. Numéro de lot. Número de lote. Chargen-Nummer. Numero di lotto. Código do lote. Partnummer. Lot nummer. Satsnummer.
	Catalogue number. Référence catalogue. Número de catálogo. Katalog-Nummer. Numero di catalogo. Número catalogo. Katalognummer.
	Use-by date. Date de péremption. Fecha de caducidad. Verfallsdatum. . La data di scadenza. Prazo de validade. Udløbsdato. Utløpsdato. Använd före.
	Temperature limit. Seuil de températures. Rango de temperature. Temperaturbereich. Limitazioni di temperatura. Limite de temperatura. Opbevaringstemperatur. Oppbevares ved. Förvaringstemperatur.
	Biological risk. Risque biologique. Riesgo biológico. Biologische Gefährdung. Rishio biologico. Risco biológico. Biologisk risk.
	Consult instructions for use. Lire le mode d'emploi. Consulte las instrucciones de uso. Gebrauchsanweisung beachten. Leggere le istruzioni per l'uso. Consultar as instruções de utilização. Se brugsanvisning. Se bruksanvisningen. Läs instruktionsmanualen.
	In vitro diagnostic medical device. Dispositif médical de diagnostic in vitro. Producto sanitario para diagnóstico in vitro. In-vitro-Diagnostikum. Dispositivo medico-diagnostico in vitro. Dispositivos médicos para diagnóstico <i>in vitro</i> . In Vitro medisinsk diagnoseutstyr. In vitro diagnostik medicinsk utrustning.
	Warning. Attention. Atensi3n. Achtung. Attenzione. Atenç3o. Advarsel. Advarsel. Varning.
	Manufacturer. Fabricant. Fabricante. Hersteller. Produttore. Fabricante. Fabrikant. Producent. Produsent. Tillverkare.
	Contains sufficient for 96 tests. Contenu suffisant pour 96 tests. Contenido suficiente para 96 pruebas. Inhalt ausreichend für 96 Tests. Contenuto sufficiente per 96 test. Indeholder tilstrækkelig for 96 test. Inneholder tilstrekkelig for 96 test. Innehåller tillräckligt för 96 test.
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive. Conformément à la directive européenne 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro. La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnóstico in vitro. Konform mit Richtlinie 98/79/EG zu In-vitro-Diagnostika. Conformità alla direttiva 98/79/CE relativa ai dispositivi medico-diagnostici in vitro. Medicinsk utstyr til in vitro-diagnostik, i overensstemmelse med Europa-Parlamentets og Rådets direktiv 98/79/EF. Medisinsk utstyr i samsvar med EU in vitro diagnostic directive 98/79/EF. Överensstämmer med direktiv 98/79/EG för medicintekniska produkter.

<b>Ag</b>	Antigen. Antigène. Antigeno, Antigene. L'antigene.
<b>DIL</b>	Diluent. Diluant. Diluyente. Probenverdünnungspuffer. Il diluente. Spädningsbuffert.
<b>CONJ</b>	Conjugate. Conjugué. Conjugado. Conjugato. Konjugat.
<b>BUF</b>   <b>WASH</b>   <b>30X</b>	Wash solution 30x conc. Solution lavage conc. 30x. Solución de lavado conc. 30x. Waschpuffer 30x konc. Soluzione di lavaggio 30x conc. Vaskebuffert/Tvättbuffert 30x konc.
<b>SUBS</b>   <b>pNPP</b>	Substrate pNPP. pNPP Substrat. Sustrato pNPP. Substrato pNPP. Substrat pNPP.
<b>CONTROL</b>   <b>X</b>	Control. Kontrolle. Contrôle. Controllo. Kontroll.

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- 1. van der Woude FJ, Rasmussen N, Lobatto S, Wiik A, Permin H, Van Es L et al.** Autoantibodies against neutrophils and monocytes: Tool for diagnosis and marker of disease activity in Wegener's granulomatosis. *Lancet* 1985, i 425-429.
- 2. Wiik A, Rasmussen N, Wieslander J.** Methods to detect autoantibodies to neutrophilic granulocytes, in Van Venrooij WJ, Maini RN (eds): *Manual of biological markers of disease*. Kluwer Academic Publishers 1993, A9, 1-14.
- 3. Rasmussen N, Sjölin C, Isaksson B, Bygren P, Wieslander J.** An ELISA for the detection of anti-neutrophil cytoplasm antibodies (ANCA). *J Immunol Methods* 1990, 127, 139-145.
- 4. Segelmark M, et al.** How and why should we detect ANCA ? *Clin Exp Rheumatol* 2000, 18, 629-635.
- 5. Jenne DE, Tschopp J, Lüdemann J, Utecht B, Gross WL.** Wegeners autoantigen decoded. *Nature* 1990, 346, 520.
- 6. Wieslander J, Wiik A.** ANCA antigens: Proteinase 3, in van Venrooij W, Maini R (eds): *Manual of Biological Markers of Disease*. Kluwer Academic Publishers 1994, 1-9.
- 7. Falk RJ, Jennette JC.** Anti-neutrophil cytoplasmic autoantibodies with specificity for myeloperoxidase in patients with systemic vasculitis and idiopathic necrotizing and crescentic glomerulonephritis. *N Engl J Med* 1988, 318, 1651-1657.
- 8. Wieslander J, Wiik A.** ANCA antigens: Myeloperoxidase, in van Venrooij W, Maini R (eds): *Manual of Biological Markers of Disease*. Kluwer Academic Publishers 1994, 1-9.
- 9. Kamech L, Harper L, Savage C.** ANCA - Positive Vasculitis. I. *Am Soc Nephrol* 2002 43 1953 - 1960.
- 10. Segelmark M, Elzouki AN, Wieslander J, Eriksson S.** Heterozygosity for the alpha 1-antitrypsin PiZ gene affects the outcome of PR3-ANCA positive vasculitis. *Kidney Int* 1995, 48, 844-850.
- 11. Boomsma MM, Stegeman CA, van der Leij MJ, et al.** Prediction of relapses in Wegener's granulomatosis by measurement of ANCA levels: a prospective study. *Arthritis Rheum*, 2000, 43 2025 - 2033.
- 12. Savige J, Davies D, Falk RJ, Jennette JC, and Wiik A.** ANCA and associated diseases: A review of the clinical and laboratory features. *Kidney Int* 2000, 57, 846-862.
- 13. Savige J, Gills D, Benson E et al.** International Consensus Statement on Testing and Reporting of ANCA. *Am J Clin Pathol* 1999, 111, 507-513.
- 14. Won H, et al.** Serial ANCA titers. Useful tool for the prevention of relapses in ANCA associated vasculitis. *Kidney Int*, 2003, 63, 1079-1085.

**Explanation of symbols. L'explication de symboles. La explicación de símbolos. Erklärung der Symbole. La spiegazione di simboli. Forklaring til symboler. Förklaringar till symboler.**

	Batch code. Numéro de lot. Número de lote. Chargen-Nummer. Numero di lotto. Código do lote. Partnummer. Lot nummer. Satsnummer.
	Catalogue number. Référence catalogue. Número de catálogo. Katalog-Nummer. Numero di catalogo. Número catalogo. Katalognummer.
	Use-by date. Date de péremption. Fecha de caducidad. Verfallsdatum. . La data di scadenza. Prazo de validade. Udløbsdato. Utløpsdato. Använd före.
	Temperature limit. Seuil de températures. Rango de temperature. Temperaturbereich. Limitazioni di temperatura. Limite de temperatura. Opbevaringstemperatur. Oppbevares ved. Förvaringstemperatur.
	Biological risk. Risque biologique. Riesgo biológico. Biologische Gefährdung. Rishio biologico. Risco biológico. Biologisk risk.
	Consult instructions for use. Lire le mode d'emploi. Consulte las instrucciones de uso. Gebrauchsanweisung beachten. Leggere le istruzioni per l'uso. Consultar as instruções de utilização. Se brugsanvisning. Se bruksanvisningen. Läs instruktionsmanualen.
	In vitro diagnostic medical device. Dispositif médical de diagnostic in vitro. Producto sanitario para diagnóstico in vitro. In-vitro-Diagnostikum. Dispositivo medico-diagnostico in vitro. Dispositivos médicos para diagnóstico <i>in vitro</i> . In Vitro medisinsk diagnoseutstyr. In vitro diagnostik medicinsk utrustning.
	Warning. Attention. Atensi3n. Achtung. Attenzione. Atenç3o. Advarsel. Advarsel. Varning.
	Manufacturer. Fabricant. Fabricante. Hersteller. Produttore. Fabricante. Fabrikant. Producent. Produsent. Tillverkare.
	Contains sufficient for 96 tests. Contenu suffisant pour 96 tests. Contenido suficiente para 96 pruebas. Inhalt ausreichend für 96 Tests. Contenuto sufficiente per 96 test. Inneholder tilstrækkelig for 96 test. Inneholder tilstrekkelig for 96 test. Innehåller tillräckligt för 96 test.
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive. Conformément à la directive européenne 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro. La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnóstico in vitro. Konform mit Richtlinie 98/79/EG zu In-vitro-Diagnostika. Conformità alla direttiva 98/79/CE relativa ai dispositivi medico-diagnostici in vitro. Medicinsk utstyr til in vitro-diagnostik, i overensstemmelse med Europa-Parlamentets og Rådets direktiv 98/79/EF. Medisinsk utstyr i samsvar med EU in vitro diagnostic directive 98/79/EF. Överensstämmer med direktiv 98/79/EG för medicintekniska produkter.

<b>Ag</b>	Antigen. Antigène. Antigeno, Antigene. L'antigene.
<b>DIL</b>	Diluent. Diluant. Diluyente. Probenverdünnungspuffer. Il diluente. Spädningsbuffert.
<b>CONJ</b>	Conjugate. Conjugué. Conjugado. Conjugato. Konjugat.
<b>BUF</b>   <b>WASH</b>   <b>30X</b>	Wash solution 30x conc. Solution lavage conc. 30x. Solución de lavado conc. 30x. Waschpuffer 30x konc. Soluzione di lavaggio 30x conc. Vaskebuffert/Tvättbuffert 30x konc.
<b>SUBS</b>   <b>pNPP</b>	Substrate pNPP. pNPP Substrat. Sustrato pNPP. Substrato pNPP. Substrat pNPP.
<b>CONTROL</b>   <b>X</b>	Control. Kontrolle. Contrôle. Controllo. Kontroll.

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**EURO DIAGNOSTICA AB**