Instruction

WIESLAB[®] ANCA panel kit

Qualitative specificity assay

Enzyme immunoassay for detection of autoantibodies
against six ANCA antigens
(Azurocidin, BPI, Cathepsin G, Elastase, Lactoferrin and Lysozyme)

Microtitration strips (12x8) 96 wells
Store the kit at +2-8° C
For in vitro diagnostic use only

Document No. E-23-0168-07
May, 2015

WIESLAB[®] ANCA panel kit Qualitative specificity assay

English:	page	2
Français:	page	11
Español:	página	18
Deutsch:	Seite	24
Italiano:	pagina	31
Dansk:	side.....	38
Norsk:	side.....	44
Svenska:	sida.....	50

REF

PAN 106

IVD

INTENDED USE

The Wieslab[®] ANCA panel test kit is an enzyme-linked immunosorbent assay (ELISA) for qualitative detection of IgG antibodies to azurocidin, BPI, cathepsin G, elastase, lactoferrin, and lysozyme in human sera. The assay is used to detect antibodies in a single serum specimen.

The test is useful as an aid in screening for specificity of ANCA IIF positive samples which are negative for PR3-ANCA and MPO-ANCA, or for detection of ANCA in inflammatory bowel diseases, liver diseases etc. One patient's sample can be tested for all six antigens on twelve occasions.

The analysis should be performed by trained laboratory professionals.

FOR IN VITRO DIAGNOSTIC USE.

Summary and explanation

ANCAs (anti-neutrophil cytoplasmic antibodies) are a family of autoantibodies related to vasculitis and inflammatory disorders. Since 1985, when c-ANCA was shown to be related to Wegener's granulomatosis (1), interest in ANCAs has steadily increased, and today these antibodies are considered to be major diagnostic tools for the investigation of systemic vasculitis.

The first method to detect ANCA was indirect immunofluorescence (IIF) performed on ethanol fixed granulocytes (2). This method yields two patterns, a cytoplasmic staining of the granulocyte denoting the presence of c-ANCAs, and a perinuclear staining denoting the presence of p-ANCAs. IIF was followed by ELISAs using the purified proteins (3, 4).

In systemic vasculitis the two most important antigens are proteinase 3 and myeloperoxidase. In inflammatory conditions like rheumatoid arthritis, autoimmune hepatitis or inflammatory bowel disease however, many other antigens have been proposed to be involved. The most important being azurocidin, BPI, (BPI=Bactericidal permeability increasing protein), cathepsin G, elastase, lactoferrin and lysozyme.

Principle of the Wieslab[®] ANCA screening assay

The wells of a microtitre strips are coated with purified ANCA antigens. During the first incubation, specific antibodies in diluted serum, will bind to the antigen coating.

The wells are then washed to remove unbound antibodies and other components.

A conjugate of alkaline phosphatase-labelled antibodies to human IgG binds to the antibodies in the wells in this second incubation.

After a further washing step, detection of specific antibodies is obtained by incubation with substrate solution. The amount of bound antibodies correlates to the colour intensity and is measured in terms of absorbance (optical density (OD)).

Warnings and precautions

- For in vitro diagnostic use.
- The human serum components used in the preparation of the controls and calibrators in the kit have been tested for the presence of antibodies to human immunodeficiency virus 1 & 2 (HIV 1&2), hepatitis C (HCV) as well as hepatitis B surface antigen by FDA approved methods and found negative. Because no test methods can offer complete assurance that HIV, HCV, hepatitis B virus, or other infectious agents are absent, specimens and human-based reagents should be handled as if capable of transmitting infectious agents.
- The Centers for Disease Control and Prevention and National Institutes of Health recommended that potentially infectious agents be handled at the Biosafety Level 2.
- All solutions contain ProClin 300 as a preservative. Never pipette by mouth or allow reagents or patient sample to come into contact with skin. Reagents containing ProClin may be irritating. Avoid contact with skin and eyes. In case of contact, flush with plenty of water.
- The concentrations of ANCA in a given specimen determined with assays from different manufacturers can vary due to differences in assay methods and reagent specificity.

Material safety data sheets for all hazardous components contained in this kit are available on request from Euro Diagnostica.

Warning

BUF	WASH	30X
DIL		

CONJ	
SUBS	pNPP

Contains ProClin 300:
 Reaction mass of: 5-chloro-2-methyl-4-isothiazolin-3-one [EC no. 247-500-7] and 2-methyl-4-isothiazolin-3-one [EC no. 220-239-6] (3:1)

- H317: May cause an allergic skin reaction.
- P264: Wash hands thoroughly after handling.
- P280: Wear protective gloves/protective clothing/eye protection/face protection.
- P302+352: IF ON SKIN: Wash with plenty of soap and water.
- P333+313: If skin irritation or rash occurs: Get medical advice/attention.

Specimen collection

The ANCA panel assay is for use with serum. Handle as if capable of transmitting agents. Avoid using sera which are icteric, lipemic and hemolyzed. Heat-inactivated sera can yield unspecific reactivities and should not be used. Store serum between 2-8° C if testing will take place within five days. If specimens are to be kept for longer periods, store at -20° C or colder. Do not use a frost-free freezer because it may allow the specimens to go through freeze-thaw cycles and degrade antibody. Samples that are improperly stored or are subjected to multiple freeze-thaw cycles may yield spurious results. The NCCLS provides recommendations for storing blood specimens, (Approved Standard-Procedures for the Handling and Processing of Blood Specimens, H18A, 1990).

Kit components and storage of reagents

- One frame with strips (12x8) coated with ANCA antigens and one lid sealed in a foil pack .
- 13 mL conjugate containing alkaline phosphatase-labelled antibodies to human IgG (blue colour).
- 32 mL Diluent (Dil) containing PBS (red colour).
- 13 mL Substrate pNPP.
- 30 mL wash solution 30x concentrated.

All reagents in the kit are ready for use except wash solution and should be stored at + 2-8° C. Remove only the number of strips needed for testing, resealing the aluminium package carefully.

Materials or equipment required but not provided

- Microplate reader with filter 405 nm.
- Precision pipettes with disposable tips.
- Washer for strips, absorbent tissue, tubes and a timer.

PROCEDURE

All solutions should be used at room temperature. One strip should be used for one patient sample. Use the enclosed protocol for ANCA panel when performing the assay.

Preparation of washing solution

In case salt crystals are observed in the vial with concentrated wash solution, place the vial at 37° C water bath until the crystals have dissolved before dilution of wash solution. Dilute 30 mL of the 30x concentrated wash solution in 870 mL distilled water. When stored at 2-8° C, the diluted wash solution is stable until the date of expiration of the kit.

Dilution of serum and incubation

Dilute the patient's serum 1/80 with diluent (790 µl diluent + 10 µl serum).

Pipet 100 µl of diluted serum in well H and 100 µl diluent in well H according to the schema below.

Column 1	Antigen coat
100 µl diluted serum	← HSA as blank
100 µl diluted serum	← Azurocidin
100 µl diluted serum	← BPI
100 µl diluted serum	← Cathepsin G
100 µl diluted serum	← Elastase
100 µl diluted serum	← Lactoferrin
100 µl diluted serum	← Lysozyme
100 µl diluent	← Human IgG (Assay Control)

Put on the lid and incubate for 30 minutes at room temperature (+20-25° C).

After serum incubation

Wash 3 times with washing solution, filling and emptying the wells each time; after the last wash, empty the wells by tapping the strip on an absorbent tissue.

Adding conjugate

Add 100 µl conjugate to each well. Put on the lid and incubate for 30 minutes in room temperature (+20-25° C).

After conjugate incubation

Wash 3 times with washing solution, filling and emptying the wells each time; after the last wash, empty the wells by tapping the strip on an absorbent tissue.

Adding substrate solution

Add 100 µl substrate pNPP solution to each well. Put on the lid and incubate for 60 minutes in room temperature (+20-25° C).

Read the absorbance at 405 nm on a microplate reader.

Calculations

The patient sample is negative for azurocidin, elastase, laktoferrin and lysozyme if the OD value is < 0,4.

The patient sample is negative for BPI and cathepsin G if the OD value is <0,6

If the patient OD value is > 0,4 on azurocidin or elastase or laktoferrin or lysozyme or > 0,6 on BPI or cathepsin G an optical density (OD) ratio for the antigen is calculated as follows:

$$\text{OD ratio} = \frac{\text{OD of patient sample for the antigen}}{\text{OD of patient sample blank (well A)}}$$

The patient sample is **negative** if the OD ratio is < 3,0 and **positive** if the OD ratio is > 4,0.

3.0 – 4,0 is **equivocal**. Retest, if still equivocal retest by an alternative method or test a new sample.

Quality Control

The OD for the patient sample in the blank well (A) should be $< 0,35$

The OD for the assay control (well H) should be $> 1,0$

If the OD value for the patient sample is > 0.35 in the blank well (A), the sample contains unspecific reactivities, which has to be checked for.

If the assay control value is not within their respective ranges, the test should be considered invalid and the test should be repeated.

The patient sample may contain unspecific reactivities if it is positive on more than one antigen.

Interpretation of results

A sample with a OD ratio of : $< 3,0 = \text{Negative}$ $> 4,0 = \text{Positive}$

$3.0 - 4,0 = \text{Equivocal}$; Retest, if still equivocal retest by an alternative method or test a new sample.

Limitations

The individual patient's OD ratio can not be used as a measure of disease severity, as antibodies from different patients may differ from each other in affinity. Thus, it is difficult to obtain an absolute standardisation of results.

The test should not be relied upon as the sole basis of decisions on clinical therapy, but should be used in combination with clinical symptoms and the results of other available tests.

Sera from patients with other autoimmune diseases and from normal individuals may contain potentially cross-reactive autoantibodies. Some individuals may be positive, with little or no evidence of clinical disease. On the other hand, some patients with active disease may have undetectable levels of these antibodies.

Immunosuppressive therapy should not be started on basis of a positive ANCA result. Initiation or changes in treatment should not be based on changes in ANCA concentration alone, but rather on careful clinical observation.

Expected results

$< 5\%$ ANCA are found in normal healthy individuals. Normally around 1% of samples sent to a laboratory contain unspecific reactivities, which have to be checked for.

Performance characteristics

Table 1. Clinical sensitivity and specificity. A total of 258 frozen retrospective sera with clinical characterisation were assayed. The following table summarises the results. (Equivocal samples are not included in the calculations).

Control and Disease groups	Total number	AZ pos	BPI pos	CH pos	EL pos	LF pos	LY pos
Blood donors:	90	0	4	4	0	0	0
Syst. Vasculitis:	48	2	2	1	0	0	0
Cystic fibrosis:	31	0	17	1	0	0	0
Ulcerative colitis:	46	2	1	1	0	0	1
Coeliac disease:	46	0	0	2	0	0	1
SCDA:	42	0	1	10	0	0	0
Crohns syndrome:	45	0	1	1	0	0	2

IBD = Inflammatory bowel diseases, **SCDA** = serum containing different autoantibodies (ANA, DNA, SSA etc.)

Table 2. Batch to batch-variation was determined by testing one sample in three -five different batches.

	AZ	BPI	CH	EL	LF	LY
Batch #	5	4	4	4	3	5
Mean OD	2,1	1,2	1,6	2,3	1,6	1,8
SD	0,3	0,15	0,35	0,25	0,03	0,49
CV %	14	12	22	11	2	13

Table 3. Inter-assay precision was determined by testing two samples in six different runs.

	AZ		BPI		CH		EL		LF		LY	
	1	2	1	2	1	2	1	2	1	2	1	2
Sample	1,5	2,6	1,2	1,3	1,6	1,3	1,8	2,1	1,6	1,1	1,7	0,7
OD	1,5	2,6	1,2	1,3	1,6	1,3	1,8	2,1	1,6	1,1	1,7	0,7
SD	0,12	0,05	0,13	0,1	0,05	0,09	0,14	0,11	0,13	0,08	0,06	0,09
CV %	8	2	12	7	3	6	8	5	8	7	4	14

Table 4. Intra-assay precision was determined by testing one sample in 12 wells.

	AZ	BPI	CH	EL	LF	LY
Sample	1	1	1	1	1	1
Mean OD	1,8	1,1	1,9	2,0	1,8	1,4
SD	0,08	0,03	0,05	0,14	0,08	0,04
CV%	5	3	2	7	5	3

Troubleshooting

Problem:	Possible causes:	Solution:
Control values out of range.	<ol style="list-style-type: none"> 1. Incorrect temperature, timing or pipetting; reagents not mixed. 2. Cross contamination of controls. 3. Improper dilution. 4. Optical pathway not clean. 	<ol style="list-style-type: none"> 1. Check that the time and temperature was correct. See "Poor precision" below. Repeat test. 2. Pipette carefully. 3. Repeat test. 4. Check for dirt or air bubbles in the wells. Wipe bottom and reread.
All test results negative.	<ol style="list-style-type: none"> 1. One or more reagents not added, or added in wrong sequence. 2. Antigen coated plate inactive. 	<ol style="list-style-type: none"> 1. Recheck procedure. Check for unused reagents. Repeat test. 2. Check for obvious moisture in unused wells. Wipe bottom and reread.
All test results yellow.	<ol style="list-style-type: none"> 1. Contaminated buffers or reagents. 2. Washing solution contaminated. 3. Improper dilution of serum. 	<ol style="list-style-type: none"> 1. Check all solutions for turbidity. 2. Use clean container. Check quality of water used to prepare solution. 3. Repeat test.
Poor precision.	<ol style="list-style-type: none"> 1. Pipette delivery CV greater than 5%. 2. Serum or reagents not mixed sufficiently or not equilibrated to room temperature. 3. Reagent addition taking too long; inconsistency in timing intervals. 4. Optical pathway not clean. 5. Washing not consistent; trapped bubbles; washing solution left in the wells. 6. Improper pipetting. 	<ol style="list-style-type: none"> 1. Check calibration of pipette. Use reproducible technique. 2. Mix all reagents gently but thoroughly and equilibrate to room temperature. 3. Develop consistent uniform technique and use multi-tip device or auto dispenser to decrease time. 4. Check for air bubbles in the wells. Wipe bottom and reread. 5. Check that all wells are filled and aspirated uniformly. Dispense liquid above level of reagent in well. After the last wash, empty the wells by tapping the strip on an absorbent tissue. 6. Avoid air bubbles in pipette tips.

REFERENCES. RÉFÉRENCES. REFERENCIAS. LITERATUR. BIBLIOGRAFIA. LITTERATUR. REFERANSER. REFERENSER

Bygren, P., Rasmussen, N., Isaksson, B., and Wieslander, J. (1992). Anti-neutrophil cytoplasm antibodies, anti-GBM antibodies and anti-dsDNA antibodies in glomerulonephritis. *Eur. J. Clin. Invest.* 22, 783-792.

Cohen Tervaert, J.W., Mulder, L., Stegeman, C., Elema, J., Huitema, M., The, H., and Kallenberg, C. (1993). Occurrence of autoantibodies to human leucocyte elastase in Wegener's granulomatosis and other inflammatory disorders. *Ann. Rheum. Dis.* 52, 115-120.

Coremans, I.E., Hagen, E.C., van der Voort, E.A., van der Woude, F.J., Daha, M.R., and Breedveld, F.C. (1993). Autoantibodies to neutrophil cytoplasmic enzymes in Felty's syndrome. *Clin. Exp Rheumatol.* 11, 255-262.

Halbwachs Mecaelli, L., Nusbaum, P., Noel, L.H., Reumaux, D., Erlinger, S., Grunfeld, J.P., and Lesavre, P. (1992). Antineutrophil cytoplasmic antibodies (ANCA) directed against cathepsin G in ulcerative colitis, Crohn's disease and primary sclerosing cholangitis. *Clin. Exp. Immunol.* 90, 79-84.

Lesavre, P., Nusbaum, P., and Halbwachs Mecaelli, L. (1993). Methods of detection of anti-cathepsin G autoantibodies in human. *Adv. Exp. Med. Biol.* 336, 257-261.

Peen, E., et al. Anti-lactoferrin antibodies and other types of ANCA in ulcerative colitis, primary sclerosing cholangitis, and Crohn's disease. *Gut.* 34, 56-62.

Skogh, T. and Peen, E. (1993). Lactoferrin, anti-lactoferrin antibodies and inflammatory disease. *Adv. Exp. Med. Biol.* 336, 533-538.

Zhao MH, Lockwood CM. Azurocidin is a novel antigen for anti-neutrophil cytoplasmic autoantibodies (ANCA) in systemic vasculitis. *Clin Exp Immunol* 1996, 103, 397-402.

Zhao MH, Jones S.J., Lockwood C.M. Bactericidal/permeability-increasing protein (BPI) is an important antigen for anti-neutrophil cytoplasmic autoantibodies (ANCA) in vasculitis. *Clin Exp Immunol* 1995; 99:49-56.

Explanation of symbols. L'explication de symboles. La explicación de símbolos. Erklärung der Symbole. La spiegazione di simboli. Forklaring til symboler. Symbolforklaring. Förklaringar till symboler.

	Batch code. Numéro de lot. Número de lote. Chargen-Nummer. Numero di lotto. Código do lote. Partnummer. Lot nummer. Satsnummer.
	Catalogue number. Référence catalogue. Número de catálogo. Katalog-Nummer. Numero di catalogo. Número catalogo. Katalognummer.
	Use-by date. Date de péremption. Fecha de caducidad. Verfallsdatum. . La data di scadenza. Prazo de validade. Udløbsdato. Utløpsdato. Använd före.
	Temperature limit. Seuils de températures. Rango de temperature. Temperaturbereich. Limitazioni di temperatura. Limite de temperatura. Opbevaringstemperatur. Oppbevares ved. Förvaringstemperatur.
	Biological risk. Risque biologique. Riesgo biológico. Biologische Gefährdung. Rishio biologico. Risco biológico. Biologisk risk.
	Consult instructions for use. Lire le mode d'emploi. Consulte las instrucciones de uso. Gebrauchsanweisung beachten. Leggere le istruzioni per l'uso. Consultar as instruções de utilização. Se brugsanvisning. Se bruksanvisningen. Läs instruktionsmanualen.
	In vitro diagnostic medical device. Dispositif médical de diagnostic in vitro. Producto sanitario para diagnóstico in vitro. In-vitro-Diagnostikum. Dispositivo medico-diagnostico in vitro. Dispositivos médicos para diagnóstico <i>in vitro</i> . In Vitro medisinsk diagnoseutstyr. In vitro diagnostik medicinsk utrustning.
	Warning. Attention. Atención. Achtung. Attenzione. Atenção. Advarsel. Advarsel. Varning.
	Manufacturer. Fabricant. Fabricante. Hersteller. Produttore. Fabricante. Fabrikant. Producent. Produsent. Tillverkare.
	Contains sufficient for 96 tests. Contenu suffisant pour 96 tests. Contenido suficiente para 96 pruebas. Inhalt ausreichend für 96 Tests. Contenuto sufficiente per 96 test. Indeholder tilstrækkelig for 96 test. Inneholder tilstrekkelig for 96 test. Innehåller tillräckligt för 96 test.
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive. Conformément à la directive européenne 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro. La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnóstico in vitro. Konform mit Richtlinie 98/79/EG zu In-vitro-Diagnostika. Conformità alla direttiva 98/79/CE relativa ai dispositivi medico-diagnostici in vitro. Medicinsk utstyr til in vitro-diagnostik, i overensstemmelse med Europa-Parlamentets og Rådets direktiv 98/79/EF. Medisinsk utstyr i samsvar med EU in vitro diagnostic directive 98/79/EF. Överensstämmer med direktiv 98/79/EG för medicintekniska produkter.

<table border="1"> <tr> <td data-bbox="395 197 491 257">Ag</td> </tr> </table>	Ag	Antigen. Antigène, Antigeno. Antigene. L'antigene. Antigen. Antigen. Antigen.		
Ag				
<table border="1"> <tr> <td data-bbox="395 300 491 360">DIL</td> </tr> </table>	DIL	Diluent. Diluant. Diluyente. Probenverdünnungspuffer. Il diluente. Fortynning. Spädningsbuffert.		
DIL				
<table border="1"> <tr> <td data-bbox="355 403 491 463">CONJ</td> </tr> </table>	CONJ	Conjugate. Conjugué. Conjugado. Conjugato. Konjugat. Konjugat. Konjugat.		
CONJ				
<table border="1"> <tr> <td data-bbox="124 528 256 589">BUF</td> <td data-bbox="256 528 405 589">WASH</td> <td data-bbox="405 528 507 589">30X</td> </tr> </table>	BUF	WASH	30X	Wash solution 30x conc. Solution lavage conc. 30x. Solución de lavado conc. 30x. Waschpuffer 30x konc. Soluzione di lavaggio 30x conc. Vaskebuffert 30x konc. Vaskeløsning 30x kons. Tvättbuffert 30x konc.
BUF	WASH	30X		
<table border="1"> <tr> <td data-bbox="220 640 352 701">SUBS</td> <td data-bbox="352 640 491 701">pNPP</td> </tr> </table>	SUBS	pNPP	Substrate pNPP. pNPP Substrat. Sustrato pNPP. Substrato pNPP. Substrat pNPP. Substrat pNPP.	
SUBS	pNPP			

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- Bygren, P., Rasmussen, N., Isaksson, B., and Wieslander, J. (1992).** Anti-neutrophil cytoplasm antibodies, anti-GBM antibodies and anti-dsDNA antibodies in glomerulonephritis. *Eur. J. Clin. Invest.* 22, 783-792.
- Cohen Tervaert, J.W., Mulder, L., Stegeman, C., Elema, J., Huitema, M., The, H., and Kallenberg, C. (1993).** Occurrence of autoantibodies to human leucocyte elastase in Wegener's granulomatosis and other inflammatory disorders. *Ann. Rheum. Dis.* 52, 115-120.
- Coremans, I.E., Hagen, E.C., van der Voort, E.A., van der Woude, F.J., Daha, M.R., and Breedveld, F.C. (1993).** Autoantibodies to neutrophil cytoplasmic enzymes in Felty's syndrome. *Clin. Exp Rheumatol.* 11, 255-262.
- Halbwachs Mecaelli, L., Nusbaum, P., Noel, L.H., Reumaux, D., Erlinger, S., Grunfeld, J.P., and Lesavre, P. (1992).** Antineutrophil cytoplasmic antibodies (ANCA) directed against cathepsin G in ulcerative colitis, Crohn's disease and primary sclerosing cholangitis. *Clin. Exp. Immunol.* 90, 79-84.
- Lesavre, P., Nusbaum, P., and Halbwachs Mecaelli, L. (1993).** Methods of detection of anti-cathepsin G autoantibodies in human. *Adv. Exp. Med. Biol.* 336, 257-261.
- Peen, E., et al.** Anti-lactoferrin antibodies and other types of ANCA in ulcerative colitis, primary sclerosing cholangitis, and Crohn's disease. *Gut.* 34, 56-62.
- Skogh, T. and Peen, E. (1993).** Lactoferrin, anti-lactoferrin antibodies and inflammatory disease. *Adv. Exp. Med. Biol.* 336, 533-538.
- Zhao MH, Lockwood CM.** Azurocidin is a novel antigen for anti-neutrophil cytoplasmic autoantibodies (ANCA) in systemic vasculitis. *Clin Exp Immunol* 1996, 103, 397-402.
- Zhao MH, Jones S.J., Lockwood C.M.** Bactericidal/permeability-increasing protein (BPI) is an important antigen for anti-neutrophil cytoplasmic autoantibodies (ANCA) in vasculitis. *Clin Exp Immunol* 1995; 99:49-56.

Explanation of symbols. L'explication de symboles. La explicación de símbolos. Erklärung der Symbole. La spiegazione di simboli. Forklaring til symboler. Symbolforklaring. Förklaringar till symboler.

	Batch code. Numéro de lot. Número de lote. Chargen-Nummer. Numero di lotto. Código do lote. Partnummer. Lot nummer. Satsnummer.
	Catalogue number. Référence catalogue. Número de catálogo. Katalog-Nummer. Numero di catalogo. Número catalogo. Katalognummer.
	Use-by date. Date de péremption. Fecha de caducidad. Verfallsdatum. . La data di scadenza. Prazo de validade. Udløbsdato. Utløpsdato. Använd före.
	Temperature limit. Seuils de températures. Rango de temperature. Temperaturbereich. Limitazioni di temperatura. Limite de temperatura. Opbevaringstemperatur. Oppbevares ved. Förvaringstemperatur.
	Biological risk. Risque biologique. Riesgo biológico. Biologische Gefährdung. Rishio biologico. Risco biológico. Biologisk risk.
	Consult instructions for use. Lire le mode d'emploi. Consulte las instrucciones de uso. Gebrauchsanweisung beachten. Leggere le istruzioni per l'uso. Consultar as instruções de utilização. Se brugsanvisning. Se bruksanvisningen. Läs instruktionsmanualen.
	In vitro diagnostic medical device. Dispositif médical de diagnostic in vitro. Producto sanitario para diagnóstico in vitro. In-vitro-Diagnostikum. Dispositivo medico-diagnostico in vitro. Dispositivos médicos para diagnóstico <i>in vitro</i> . In Vitro medisinsk diagnoseutstyr. In vitro diagnostik medicinsk utrustning.
	Warning. Attention. Atención. Achtung. Attenzione. Atenção. Advarsel. Advarsel. Varning.
	Manufacturer. Fabricant. Fabricante. Hersteller. Produttore. Fabricante. Fabrikant. Producent. Produsent. Tillverkare.
	Contains sufficient for 96 tests. Contenu suffisant pour 96 tests. Contenido suficiente para 96 pruebas. Inhalt ausreichend für 96 Tests. Contenuto sufficiente per 96 test. Indeholder tilstrækkelig for 96 test. Inneholder tilstrekkelig for 96 test. Innehåller tillräckligt för 96 test.
	Conformity to 98/79/EC on In Vitro Diagnostic Medical Device Directive. Conformément à la directive européenne 98/79/CE relative aux dispositifs médicaux de diagnostic in vitro. La conformidad con la Directiva 98/79/CE sobre productos sanitarios para diagnóstico in vitro. Konform mit Richtlinie 98/79/EG zu In-vitro-Diagnostika. Conformità alla direttiva 98/79/CE relativa ai dispositivi medico-diagnostici in vitro. Medicinsk utstyr til in vitro-diagnostik, i overensstemmelse med Europa-Parlamentets og Rådets direktiv 98/79/EF. Medisinsk utstyr i samsvar med EU in vitro diagnostic directive 98/79/EF. Överensstämmer med direktiv 98/79/EG för medicintekniska produkter.

<table border="1"> <tr> <td data-bbox="395 197 491 255">Ag</td> </tr> </table>	Ag	<p>Antigen. Antigène, Antigeno. Antigene. L'antigene. Antigen. Antigen. Antigen.</p>		
Ag				
<table border="1"> <tr> <td data-bbox="395 295 491 353">DIL</td> </tr> </table>	DIL	<p>Diluent. Diluant. Diluyente. Probenverdünnungspuffer. Il diluente. Fortynning. Spädningsbuffert.</p>		
DIL				
<table border="1"> <tr> <td data-bbox="355 394 491 452">CONJ</td> </tr> </table>	CONJ	<p>Conjugate. Conjugué. Conjugado. Conjugato. Konjugat. Konjugat. Konjugat.</p>		
CONJ				
<table border="1"> <tr> <td data-bbox="124 521 256 580">BUF</td> <td data-bbox="256 521 405 580">WASH</td> <td data-bbox="405 521 512 580">30X</td> </tr> </table>	BUF	WASH	30X	<p>Wash solution 30x conc. Solution lavage conc. 30x. Solución de lavado conc. 30x. Waschpuffer 30x konc. Soluzione di lavaggio 30x conc. Vaskebuffert 30x konc. Vaskeløsning 30x kons. Tvättbuffert 30x konc.</p>
BUF	WASH	30X		
<table border="1"> <tr> <td data-bbox="220 636 352 694">SUBS</td> <td data-bbox="352 636 491 694">pNPP</td> </tr> </table>	SUBS	pNPP	<p>Substrate pNPP. pNPP Substrat. Sustrato pNPP. Substrato pNPP. Substrat pNPP. Substrat pNPP.</p>	
SUBS	pNPP			

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EURO DIAGNOSTICA AB